**ROUSSELET**  
**ROBATEL**

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## Laboratory mixer-settlers

*These mixer-settlers are extremely useful in the laboratory and pilot plant for performing the studies required to adapt the theoretical principles of liquid-liquid extraction to industrial processes.*

*The following operations are performed in one stage of a mixer-settler :*

- 1/ *Intimate contact in the mixer between the multi-component feed solution, A, and the immiscible extraction solvent, B, which is selective for a given element in the feed.*
- 2/ *Gravity separation of the resulting emulsion into two phases in a settling chamber. The two phases then pass countercurrently to contiguous stages.*

*The juxtaposition of several such stages forms a mixer-settler battery.*

*The successive mixing and separation operations permit the mass transfer of the solute(s) from the feed to the solvent, producing the enriched solvent (extract), B', and the impoverished phase (raffinate), A'.*

*Since each stage is separated from contiguous stages, there is no cross-contamination when taking liquid samples. This facilitates the analysis of the equilibrium distribution of the solute(s) throughout the battery.*

### Main features

*Each stage is equipped with an interchangeable motor which drives a mixing and pumping turbine whose height can be adjusted. This turbine draws the two phases from the settlers of the adjacent stages, mixes them, and transfers this emulsion to the associated settler.*

*Hydraulic stability is ensured by the overflow weirs which regulate the inter-stage flow.*

*The overflow weirs of the heavier phase are adjustable, thus allowing the position of the interphase in each settler to be independently regulated according to the densities of the two phases.*

*The stages are machined from solid blocks of raw materials and grouped in modules. These interchangeable modules are assembled into batteries with the required number of stages.*

*Each module includes :*

- 4 stages for sizes 1 and 2 (units of 1 or 2 stages are also available in these sizes on request),
- 1 stage for sizes 3, 4, 5 and 6.

*The standard electric motors are three-phase, totally enclosed (IP 44 protection) and are very reliable since they do not contain brushes. All motors of a battery, or of a section of a battery, operate at the same speed which is adjusted by an electrical control cabinet supplying variable frequency current. The standard control cabinet is capable of regulating 8 motors simultaneously.*

*On request, other types of motors may be employed: direct current, explosion-proof (only for large mixer-settlers), or pneumatic motors.*

### Construction options

**Sizes :** Six sizes of mixer-settlers are available with many options for mixers and settlers of different volumes.

**Phase recycling :** The mixing efficiency or the phase separation may often be improved by artificially increasing the flow rate of one of the phases by recycling that phase from the settler to the mixer of the same stage.

*The internal recycling models and the universal models are equipped with recycling channels, a recycle phase selector, and a recycle flow regulator. The use of external pumps permits a more precise adjustment of the recycle flow rate.*

*The standard size 3 to 6 mixer-settlers do not have recycling channels, but they are equipped with connections for external recycling. The recycle flow rate can be adjusted by pumps, valves or orifices.*

**Sight glasses at the end of the settlers :** The internal recycling and universal models are always equipped with sight glasses. They are optional on the other models and are recommended.



**Intermediate inlets and outlets :** In addition to the internal recycling capability, the universal models have additional inlets and outlets on the mixer and settler of each stage.

*Various selectors, plugs, and connectors permit operation of the required number of stages in a given battery and allow the division of the battery into several sections for extraction, scrubbing, stripping, etc. The number of stages in each section can be adjusted according to the process requirements.*

**Covers :** Optional covers which reduce solvent evaporation may be supplied, but they are not totally vaportight. They can be made from the same material as the mixer-settler body, from plexiglass, or from glass (for sizes 1 or 2 only).

**Materials of construction :**

- polyethylene,
- polypropylene,
- polyvinylidene fluoride (for sizes 1 and 2 only),
- other materials on request (for example teflon for sizes 1 and 2).

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