



ROUSSELET ROBATEL

LABORATORY / BENCH TOP CENTRIFUGE *RA20Vx*

*For solid / liquid sedimentation
or solid / liquid centrifugal filtration*

Ideal for R & D, glass labs, chemistry labs, pharmaceutical research, and food applications...



- ❖ Compact design
- ❖ Available in stainless steel and hastelloy
- ❖ Variable speed with speed display
- ❖ Feed pipe, optional wash pipe, sightglass

Centrifuge Specifications:

Weight: 55 kg
110 V or 220 V power
Dimensions: 460 x 497 x 442 mm



Basket			
Diameter - mm	200	120	120
Useful Volume - liters	1.2	0.6	0.3
Maximum Load - kg	1.5	0.75	0.38
Maximum Speed - RPM	3000	4000	4000

ROUSSELET CENTRIFUGATION SA

Headquarters
Avenue Rhin et Danube
PAE de Marenton
07104 ANNONAY - FRANCE



+ 33 (0) 4 75 69 22 11



+ 33 (0) 4 75 67 69 80

E-mail rousselet.sa@rousselet.com

Web

www.rousselet.com

UK :
ROUSSELET ROBATEL UK Ltd
Parkside House, 17 East Parade
HARROGATE
NORTH YORKSHIRE HG1 5LF

+ 44(0) 1 423 530 093

+ 44(0) 1 423 530 120

sales@rousselet-robatel.co.uk

Germany :
ROUSSELET CENTRIFUGATION SA

Bierbacher Str. 30
D-66424 HOMBURG-WÖRSCHWEILER

+ 49(0) 6848 730580

+ 49(0) 6848 730581

rousselet@t-online.de

www.rousselet-robatel.com

USA :
ROBATEL Inc.
703 West Housatonic Street
PITTSFIELD
MA 01201

+ 1 413 499 4818

+ 1 413 499 5648

sales@robatel.com