

# FCPC® FAST CENTRIFUGAL PARTITION CHROMATOGRAPH

from LAB to INDUSTRIAL PRODUCTION

## FCPC C TECHNICAL DATA



### NEW EQUIPMENT

25 ml or 50 ml column

Very fast separation

Very low sample quantity required



DATA	25ml	50ml
Sample quantity	10 - 500 mg	10 - 1000 mg
Flow rate	1- 10 ml/min	1-10 ml/min
Separation time	5 – 20 min	
Pressure max.	80 bars / 1160 psi	
Rotation speed max.	3000 rpm	
Size / weight Power	300 x 300 x h 555 mm / 40 kg 750W – 110 or 220 VAC	
Material	Rotor: stainless steel + PTFE. Wetted parts and frame : stainless steel	

All Kromaton equipment is suitable for industrial/laboratory service, with reliable and robust operation, and is compatible with current European standards.